## 迈晨科技 M&C GENE TECHNOLOGY

PRODUCT DATASHEET

## Crystal Violet, 0.2% (w/v)

Cell Staining

产品编号: CD002 中文名称: 结晶紫溶液

产品简介: A very useful reagent for obtaining quantitative information about the relative density of cells adhering to multi-well cluster dishes. The dye in this assay, crystal violet, stains DNA. Upon solubilization, the amount of dye taken up by the monolayer can be quantitated in a spectrophotometer or plate reader. Crystal violet is an important component in Gram staining.

其它名称: Basic Violet 3, Gentian Violet, Hexamethylpararosaniline chloride, Methyl Violet 10B.

分子式: C25H30N3Cl

分子量: 407.98

CAS 编号: 548-62-9

产品类型: Ready-to-use solution

**姓** 度: Dye content in powder, ≥90%

浓 度: 0.2% (w/v)

包装规格: 100 ml

储存温度: RT

染色流程:Crystal Violet Staining can be used in many applications. The procedure provided here is for Focus Formation Assay.

- ♦ Place plates on ice.
- ♦ Wash two times with ice-cold 1X PBS.
- ♦ Fix cells with ice-cold methanol (stored at -20°C) for 10 minutes.
- Aspirate methanol from plates, move off ice and add enough 0.2% crystal violet solution to cover bottom of plate.
- ♦ Incubate at room temperature for 10 minutes.
- ♦ Pour off crystal violet solution into a beaker (can re-use).
- ♦ Carefully rinse plates in ddH20 until color no longer coming off in rinse.
- Allow plates to dry at room temperature (may take overnight).



