

Crystal Violet, 0.2% (w/v)

Cell Staining

产品编号: CD002

中文名称: 结晶紫溶液

产品简介: This is a very useful reagent for obtaining both qualitative and quantitative information on the relative density of adherent cells in a multiwell culture dish. The dye crystalline violet in this reagent stains DNA. Upon lysis, the amount of dye absorbed by a monolayer of cells can be quantified by a spectrophotometer or plate reader. Crystal violet is an important component of Gram stain.

其它名称: Basic Violet-3, Gentian Violet, Hexamethylpararosaniline chloride, Methyl Violet 10B.

分子式: C25H30N3CI

分子量: 407.98

CAS 编号: 548-62-9

产品类型: Ready-to-use solution

浓 度: 0.2% (w/v)

包装规格: 100 ml

储存温度: Room temperature (R.T.)

H₃C N CH₃

CH₃

染色流程:

Crystal Violet Staining can be used for a variety of applications. The procedure provided herein is used for Focus Formation Assay.

- 1. Crystal Violet Staining can be used in many applications. The procedure provided here is for Focus Formation Assay.
- 2. Place plates on ice.
- 3. Wash two times with ice-cold 1X PBS.
- 4. Fix cells with ice-cold methanol (stored at -20C) for 10 minutes.
- 5. Aspirate methanol from plates, move off ice and add enough 0.2% crystal violet solution to cover bottom of plate.
- 6. Incubate at room temperature for 10 minutes.
- 7. Pour off crystal violet solution into a beaker (can re-use).
- 8. Carefully rinse plates in ddH20 until color no longer coming off in rinse.
- 9. Allow plates to dry at room temperature (may take overnight).

FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC AND THERAPEUTIC PROCEDURES

MACGENE Biotechnology ● Phone: (010)8205-7786 ● (010)6237-9789

E-mail: order@macgenes.com ● Tech Support: support@macgenes.com ● URL: http://www.macgenes.com