

RNase Terminator (Spray)

Laboratory Contamination Control

Catalogue #	LS102
Description	A surface decontamination solution that destroys RNases on contact. Spray RNase Terminator Solution onto the surface to be decontaminated and rinse it off with RNase-free water. RNase Terminator can completely remove RNase contamination from glass, plastic and metal surfaces.
Properties	A low-foaming, heavy duty alkaline solution.
Formulation	Detergents and metal ion chelating agents formulated in alkali water.
Concentration	NA
Pack Size	500 ml
Stability	Store at room temperature, stable for one year from the date of shipment.
Application	Ideal solution for cleaning bench-top surfaces, pipettors, and equipment that needs to be RNase-free.
Instruction	Apply RNase Terminator Solution directly to surface to be cleaned, wipe thoroughly with paper towel, rinse with water twice and then dry with clean paper towel.
Reactivity	Skin and eye irritation/Metal surface corrosion.

FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC AND THERAPEUTIC PROCEDURES



MACGENE Biotechnology • Phone: (858) 876-8815 (USA), +86-10-8205-7786 (Asia), +381-64-26-77-772 (Europe)

E-mail: order@macgene.com • Tech Support: support@macgene.com • URL: <http://www.macgene.com>