



## Silver Staining Solution, 100X

Catalogue #: MP004  
Storage: RT  
Size: 500mL

### *Description:*

### *Notes:*

1. Do not use the plastic tray because most plastic containers have softeners (sulphur-compunds) which reduce silver nitrate to metallic silver.
2. Avoid direct contact of the gel with hands, gloves, which will leave "finger prints" on your gel after developing.
3. Prepare the working solution freshly prior to use.

### *References:*

1. Swain and Ross. Electrophoresis. 1995 Jun;16(6):948-51.
2. Blum et al. Electrophoresis. 1987; 8:93-99.

**For resaech use only**