



Iscove's Modification of DMEM (IMDM) *Cell Culture Medium*

产品编号: CM10016

产品简介: Iscove 等人的研究表明红细胞和巨噬细胞的前体细胞 (precursor cell) 能够在完全无血清的培养基中生长, 该培养基需要辅以 albumin, transferrin, lecithin 以及 selenium。IMDM 是在 DMEM 的基础上进行了改进, 加入 selenium, 氨基酸和维生素, sodium pyruvate, HEPES, 以及 potassium nitrate 取代 DMEM 中的 ferric nitrate。进一步研究表明 IMEM 支持鼠 B 淋巴细胞、骨髓中的造血组织 (hemopoietic tissue)、lipopolysaccharide 刺激的 B 细胞、T 淋巴细胞及多种杂交瘤细胞等的生长。

该产品使用 HEPES 缓冲体系, 含有 L-glutamine, 不含 alpha-thioglycerol 和 beta-mercaptoethanol。

产品类型: 无菌过滤即用型液体培养基

包装规格: 500 mL

储存条件: 2-8°C

渗透压: 285±20 mOsm

酸碱度: 7.2±0.2

参考文献:

1. Iscove, N.N and Melchers, F. (1978). Complete Replacement of Serum by Albumin, Transferrin, and Soybean Lipid in Cultures of Lipopolysaccharide-Reactive B Lymphocytes. J. Exp. Medicine. 147, 923-933.
2. Iscove, N.N., Guilbert, L.J. and Weyman, C. (1980). Complete Replacement of Serum in Primary Cultures of Erythropoietin Dependent Red Cell Precursors [CFU-E] by Albumin, Transferrin, Iron, Unsaturated Fatty Acid, Lecithin and Cholesterol. Exp. Cell Research. 126, 121-126.

FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC AND THERAPEUTIC PROCEDURES

MACGENE Biotechnology ● Phone: (010)8205-7786 ● (010)6237-9789

E-mail: order@macgenes.com ● Tech Support: support@macgenes.com ● URL: <http://www.macgene.com>



PRODUCT DATASHEET

Catalogue Number	CM15016	CM011	CM10016
Inorganic Salts	mg/L	mg/L	mg/L
CaCl ₂ (anhydrous)	165.00	165.00	165.00
KCl	330.00	330.00	330.00
KNO ₃	0.076	0.076	0.076
MgSO ₄ (anhydrous)	97.70	97.70	97.70
Na ₂ SeO ₃	0.0173	0.0173	0.0173
NaCl	4505.00	4505.00	4505.00
NaH ₂ PO ₄ •H ₂ O	125.00	125.00	125.00
NaHCO ₃	3024.00	3024.00	3024.00
Amino Acids			
L-Alanine	25.00	25.00	25.00
L-Arginine•HCl	84.00	84.00	84.00
L-Asparagine•H ₂ O	28.40	28.40	28.40
L-Aspartic acid	30.00	30.00	30.00
L-Cystine•2HCl	91.24	91.24	91.24
L-Glutamic acid	75.00	75.00	75.00
L-Glutamine			584.00
Glycine	30.00	30.00	30.00
L-Histidine•HCl•H ₂ O	42.00	42.00	42.00
L-Isoleucine	105.00	105.00	105.00
L-Leucine	105.00	105.00	105.00
L-Lysine•HCl	146.00	146.00	146.00
L-Methionine	30.00	30.00	30.00
L-Phenylalanine	66.00	66.00	66.00
L-Proline	40.00	40.00	40.00
L-Serine	42.00	42.00	42.00
L-Threonine	95.00	95.00	95.00
L-Tryptophan	16.00	16.00	16.00
L-Tyrosine•2Na•2H ₂ O	103.79	103.79	103.79
L-Valine	94.00	94.00	94.00
Vitamins			
Biotin	0.013	0.013	0.013
D-Calcium pantothenate	4.00	4.00	4.00
Choline chloride	4.00	4.00	4.00
Folic acid	4.00	4.00	4.00
i-Inositol	7.20	7.20	7.20
Nicotinamide	4.00	4.00	4.00
Pyridoxine•HCl	4.00	4.00	4.00
Riboflavin	0.40	0.40	0.40
Thiamine•HCl	4.00	4.00	4.00
Vitamin B12	0.013	0.013	0.013
Other			
D-Glucose	4500.00	4500.00	4500.00
HEPES	5958.00	5958.00	5958.00
Phenol red, Na	15.00		15.00
Sodium pyruvate	110.00	110.00	110.00

FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC AND THERAPEUTIC PROCEDURES

MACGENE Biotechnology ● Phone: (010)8205-7786 ● (010)6237-9789

E-mail: order@macgenes.com ● Tech Support: support@macgenes.com ● URL: <http://www.macgene.com>